

**KITAZATO®**

# Now you have more time to *freeze time*



## The new VT 600

*Vitrification and Thawing media for the Cryotop® Method  
Now with 12 months shelf life.*

In Kitzato we operate under the **Constant Improvement Philosophy**. The constant feedback from IVF professionals is essential to us.

**A brand that listens to its clients is a brand that improves its methods.**

We are pleased to launch the **New VT600 series**, the media that answers your storage needs with a **shelf life of 12 months**.

So, no need to rush, you can keep those vitrification and thawing procedures comming!

## Evidence of success

The Cryotop® Method media has evolved to continue offering outstanding results.

Our Method is backed up by numerous independent clinical studies with the biggest sample

### Oocytes:



n=3286  
Cobo  
2010



n=124  
Rienzi  
2010

### Blastocysts:



n=2543  
Shi  
2012



n=6861  
Liu  
2013

### Embryos:



n=2031  
Cobo  
2012



n=553  
Rodríguez-Purata  
2010



## Choose Science

Kitazato always lead vitrification with innovation. This is what makes our media special.

- **We were the first ones to replace sucrose with trehalose**, demonstrating its superiority as an osmotic agent.
- **Our media has evolved replacing the HSA (Human Serum Albumin) for HPC (Hydroxypropyl Cellulose)**, reducing the mechanical stress during thawing.
- **The incorporation of Gentamicin** prolongs the sterility of the media, guaranteeing greater safety in handling the solutions.
- **The HEPES tampon and the sodium bicarbonate act as stabilizers for the pH**, ensuring normal cellular physiological conditions durin vitrification and thawing procedures.