



SYNGA
TOUCHING LIFE GENTLY

Micromanipulation Pipettes for Embryology

For use with a micromanipulator in IVF. Sterile, MEA tested, non toxic, safe for living cells. Individually sealed pipettes in a box.

order number product

ICSI (Injection) SYNGA pipettes for IVF

ICSI Pipette with spike, very stable diameter on the whole tip to elbow length, used for sperm aspiration and intracytoplasmic injection of a single sperm into the oocyte. Blue colour coding.

SG-X30S-10 ICSI Pipette with spike, OD 6 µm, ID 5 µm, angle 30°, pack: 10x1 pipette
SG-X35S-10 ICSI Pipette with spike, OD 6 µm, ID 5 µm, angle 35°, pack: 10x1 pipette

Holding SYNGA pipettes for IVF

Holding Pipette with a rounded tip and polished opening, used to hold the oocyte or embryo when performing ICSI or other micromanipulation procedures, suitable also for blastocyst biopsy. Red colour coding.

SG-H30a-10 Holding Pipette, OD 100 µm, ID 18 µm, angle 30°, pack: 10x1 pipette
SG-H35a-10 Holding Pipette, OD 100 µm, ID 18 µm, angle 35°, pack: 10x1 pipette
SG-H30b-10 Holding Pipette for Biopsy, OD 120 µm, ID 21 µm, angle 30°, pack: 10x1 pipette
SG-H35b-10 Holding Pipette for Biopsy, OD 120 µm, ID 21 µm, angle 35°, pack: 10x1 pipette

Biopsy SYNGA pipettes for IVF

Biopsy Pipette for Blastocyst Biopsy, Day 5-6, Trophectoderm. Yellow colour coding.

SG-B30Ta-10 Blastocyst Biopsy, Day 5-6, Trophectoderm, ID 25 µm, angle 30°, pack: 10x1 pipette
SG-B35Ta-10 Blastocyst Biopsy, Day 5-6, Trophectoderm, ID 25 µm, angle 35°, pack: 10x1 pipette

Synga s.r.o. Voděradská 2552/16 251 01 Říčany u Prahy Czech Republic
VAT: CZ27923835 Tel: +420 735 613 904
info@synga.cz e-shop: www.synga.cz/eshop www.synga.cz