

Denudation and Manipulation SYNGA pipettes for IVF

CE certified, sterile, MEA tested, non toxic, safe for living cells.

Synga Flex

SG Flex - plastic pipettes

Plastic pipettes, sterile, for use with Stripper, 5x10 pipettes - 50 pipettes in a pack

Manipulation and Blastocyst transfer

SF-140-50	inner Ø 140µm, yellow
SF-180-50	inner Ø 180µm, yellow
SF-250-50	inner Ø 250µm, yellow
SF-600-50	inner Ø 600µm

Accessories

SG-007-CUT-R	SG Cutter. Storage of unopened sealed pipettes and the cutter tool to be attached to the right side wall of the IVF cabinet
SG-007-CUT-L	SG Cutter. Storage of unopened sealed pipettes and the cutter tool to be attached to the left side wall of the IVF cabinet

Samples

SG Pipettes 10 pieces

SF-140-10	inner Ø 140µm
SF-180-10	inner Ø 180µm
SF-250-10	inner Ø 250µm
SF-600-10	inner Ø 600µm