

## Denudation and Manipulation SYNGA pipettes for IVF

CE certified, sterile, MEA tested, non toxic, safe for living cells.

### I. Classic System

#### SG pipettes without bulbs (bulbs sold separately):

*SG Pipette, sterile, pack: 5x10 pipette - 50 pipettes in a pack*

##### Denudation

*Short – functional part 14 mm*

SG-130S-50	inner Ø 130µm, red
SG-140S-50	inner Ø 140µm, green
SG-145S-50	inner Ø 145µm, green
SG-150S-50	inner Ø 150µm, white
SG-160S-50	inner Ø 160µm, white

##### Denudation

*Long – functional part 29 mm*

SG-130L-50	inner Ø 130µm, red
SG-140L-50	inner Ø 140µm, green
SG-150L-50	inner Ø 150µm, white
SG-160L-50	inner Ø 160µm, white

##### Manipulation and Blastocyst transfer

SG-180-50	inner Ø 180µm, yellow
SG-250-50	inner Ø 250µm, blue

##### Oocyte Vitrification Pipettes – polished tip

SG-140P-50	inner Ø 140µm
SG-150P-50	inner Ø 150µm

### Classic Bulbs

*Bulbs - sterile, 4x5 pcs, sealed in a blister pack by 5 pieces, color may be used to indicate the diameter*

#### A. XtraSoft

SG-012-20XG	green
SG-012-20XY	yellow
SG-012-20XB	blue
SG-012-20XR	red
SG-012-20XW	white

#### B. Soft

SG-012-20SG	green
SG-012-20SY	yellow
SG-012-20SB	blue
SG-012-20SR	red
SG-012-20SW	white

**SYNGA s.r.o.**

Voděradská 2552, 251 01 Říčany u Prahy, Czech republic  
[mail@synga.cz](mailto:mail@synga.cz), [www.synga.cz](http://www.synga.cz)

## VI Bulbs

*Bulbs with reduced inner volume - sterile, 4x5 pcs, sealed in a blister pack by 5 pieces*

SG-VI-20XG	green
SG-VI-20XY	yellow
SG-VI-20XB	blue
SG-VI-20XR	red
SG-VI-20XW	white

## Accessories

SG-007-CUT-R	<b>SG Cutter.</b> Storage of unopened sealed pipettes and the cutter tool to be attached to the right side wall of the IVF cabinet
SG-007-CUT-L	<b>SG Cutter.</b> Storage of unopened sealed pipettes and the cutter tool to be attached to the left side wall of the IVF cabinet
SG-011-1	<b>SG Rack.</b> For storage of up to 6 open pipettes with bulbs in the working area in horizontal or vertical position
SG-021-3	<b>SG Pen Holder VI</b>
SG-021-4	<b>SG Pen Holder C</b>

## II. Premounted system - bulb&pipette

*SG Pipette pack: 4x5 pipette with bulb attached, blister sheets of 5 pipettes, 4 in a box, sterile single packed*

### A. Classic Bulb

SG-130S-20X-BL	inner Ø 130µm, red bulb, short pipette
SG-130L-20X-BL	inner Ø 130µm, red bulb, long pipette
SG-140S-20X-BL	inner Ø 140µm, green bulb, short pipette
SG-140L-20X-BL	inner Ø 140µm, green bulb, long pipette
SG-150S-20X-BL	inner Ø 150µm, white bulb, short pipette
SG-150L-20X-BL	inner Ø 150µm, white bulb, long pipette
SG-180-20X-BL	inner Ø 180µm, yellow bulb, short pipette
SG-250-20X-BL	inner Ø 250µm, blue bulb, short pipette

### B. VI Bulb

SG-120S-20XVI-BL	inner Ø 120µm, red bulb, short pipette
SG-120L-20XVI-BL	inner Ø 120µm, red bulb, long pipette
SG-130S-20XVI-BL	inner Ø 130µm, red bulb, short pipette
SG-130L-20XVI-BL	inner Ø 130µm, red bulb, long pipette
SG-140S-20XVI-BL	inner Ø 140µm, green bulb, short pipette
SG-140L-20XVI-BL	inner Ø 140µm, green bulb, long pipette
SG-150S-20XVI-BL	inner Ø 150µm, white bulb, short pipette
SG-150L-20XVI-BL	inner Ø 150µm, white bulb, long pipette
SG-180-20XVI-BL	inner Ø 180µm, yellow bulb, short pipette
SG-250-20XVI-BL	inner Ø 250µm, blue bulb, short pipette

## Accessories

SG-31-1	<b>SG Magazine.</b> Storage of premounted pipettes with bulbs, to be attached on the wall
---------	---

## SYNGA s.r.o.

Voděradská 2552, 251 01 Říčany u Prahy, Czech republic  
[mail@synga.cz](mailto:mail@synga.cz), [www.synga.cz](http://www.synga.cz)

### III. SAS - Synga Adhering System

#### Classic Bulbs with Adhering ball

*Classic Bulbs including metal self Adhering ball – in blister pack by 5 pieces:*

##### A. XtraSoft

SA-012-5XG	green
SA-012-5XY	yellow
SA-012-5XB	blue
SA-012-5XR	red
SA-012-5XW	white

##### B. Soft

SA-012-5SG	green
SA-012-5SY	yellow
SA-012-5SB	blue
SA-012-5SR	red
SA-012-5SW	white

#### Accessories

SA-011-1	<b>SA Magnetic Rack</b> , for storage of up to 6 open pipettes with SAS bulbs in the working area in horizontal or vertical position.
SA-021-2	<b>SA Pen Holder C</b> , Conical with magnetic dot - For manual holding with hole for push pull bulb and convenient manipulation, <b>One</b> sided, Dural, used with classical bulbs without metal ball
SA-008-R	<b>SA M dots Wall Rack</b> for Adhesive hanging up to 3 open pipettes with SAS bulbs or with a SAS Pen Holder. - <b>Right</b> side
SA-008-L	<b>SA M dots Wall Rack</b> for Adhesive hanging up to 3 open pipettes with SAS bulbs or with a SAS Pen Holder. - <b>Left</b> side
SA-003	<b>SA Base M dot Adhesive hang Rack</b> for storage of up to 3 open pipettes with SAS bulbs or pen holders

#### **IV. Samples**

**SG Pipettes and Bulbs, 10x pipette + 1x Bulb S + 1x Bulb XS + 1x Bulb VI**

SG-130S-10XS+VI	inner Ø 130µm, red bulbs, short pipette
SG-140S-10XS+VI	inner Ø 140µm, green bulbs, short pipette
SG-140L-10XS+VI	inner Ø 140µm, green bulbs, long pipette
SG-150S-10XS+VI	inner Ø 150µm, white bulbs, short pipette
SG-150L-10XS+VI	inner Ø 150µm, white bulbs, long pipette
SG-160S-10XS+VI	inner Ø 160µm, white bulbs, short pipette
SG-180-10XS+VI	inner Ø 180µm, yellow bulbs, short pipette
SG-250-10XS+VI	inner Ø 250µm, blue bulbs, short pipette

**Premounted – pipettes with bulb attached, blister of 5 pieces**

SG-130S-5XVI-BL	inner Ø 130µm, red bulb, short pipette
SG-130L-5XVI-BL	inner Ø 130µm, red bulb, long pipette
SG-140S-5XVI-BL	inner Ø 140µm, green bulb, short pipette
SG-140L-5XVI-BL	inner Ø 140µm, green bulb, long pipette
SG-150S-5XVI-BL	inner Ø 150µm, white bulb, short pipette
SG-150L-5XVI-BL	inner Ø 150µm, white bulb, long pipette
SG-180-5XVI-BL	Long narrow functional part, inner Ø 180µm, yellow bulb, short pipette
SG-250-5XVI-BL	Long narrow functional part, inner Ø 250µm, blue bulb, short pipette