



Synga Accessories

Denudation and Manipulation System for Streamlined Routines in IVF Labs

CE 2265

ISO 13485:2016



Denudation and Manipulation SYNGA pipettes for IVF

CE certified, sterile, MEA tested, non toxic, safe for living cells.

SG Accessories

SG-021-3	SG Pen Holder VI for manual holding of pipette with a VI bulb
SG-021-3-Red	SG Pen Holder VI for manual holding of pipette with a VI bulb, red
SG-021-4	SG Pen Holder C for manual holding of pipette with a classic bulb
SG-022-1	SG Ergo Pen Holder C for manual holding of pipette with a classic bulb, extra light with improved shape
SG-011-1	SG Rack for storage of up to 6 open pipettes with bulbs in the working area in horizontal or vertical position



SG Classic version

SG-007-CUT-R	SG Cutter Right for storage of unopened pipettes and the cutter tool to be attached to the right side wall of the IVF cabinet
SG-007-CUT-L	SG Cutter Left for storage of unopened pipettes and the cutter tool to be attached to the left side wall of the IVF cabinet



SA - Synga Adhering System

SA-011-1	SA Magnetic Rack for storage of up to 6 opened pipettes with SA bulbs in the working area in horizontal or vertical position.
SA-021-2	SA Pen Holder , conical with a magnetic dot, for convenient manipulation, used with classic bulbs without metal balls
SA-008-R	SA M dots Wall Rack for adhesive suspension of up to 3 opened pipettes with SA bulbs or with a SA Pen Holder. Right side
SA-008-L	SA M dots Wall Rack for adhesive suspension of up to 3 opened pipettes with SA bulbs or with a SA Pen Holder. - Left side
SA-003	SA Base M dot Adhesive hang Rack for storage of up to 3 opened pipettes with SA bulbs or pen holders



SG Clip racks

SG-003-D-G	SG Rack - Daisy Clip , with glass stand, to be placed on a desk
SG-003-D-M	SG Rack - Daisy Clip , with marble stand, to be placed on a desk
SG-003-D-A	SG Rack - Daisy Clip , with magnetic stand, to be placed on a desk
SG-003-V-G	SG Rack - V-Clip , with glass stand, to be placed on a desk
SG-003-V-M	SG Rack - V-Clip , with marble stand to be placed on a desk
SG-003-V-A	SG Rack - V-Clip , with magnetic stand to be placed on a desk
SG-008-V	SG Wall Rack - V-Clip , to be attached to a wall
SG-008-A	SG Wall Rack - Arch Clip , to be attached to a wall