



## Micromanipulation Pipettes for Embryology

For use with a micromanipulator in IVF. Sterile, MEA tested, non toxic, safe for living cells.  
Individually sealed pipettes in a box.

Revision: 2-2025

### **Synga Micromanipulation Pipettes**

*For use with a micromanipulator in IVF. Sterile, MEA, non-toxic, safe for living cells. Individually sealed pipettes in a box.*

#### **ICSI Injection SG Pipettes**

|            |   |
|------------|---|
| SG-X30S-10 | ICSI Pipette with spike, OD 5-7µm, ID 4-5 µm, angle 30°, pack: 10x1 pipette |
| SG-X35S-10 | ICSI Pipette with spike, OD 5-7µm, ID 4-5 µm, angle 35°, pack: 10x1 pipette |

#### **TESE SG Pipettes**

|            |  |
|------------|--|
| SG-X30T-10 | ICSI Pipette TESE, no spike, OD 7-9µm, ID 5-7µm, angle 30°, pack: 10x1 pipette |
| SG-X35T-10 | ICSI Pipette TESE, no spike, OD 7-9µm, ID 5-7µm, angle 35°, pack: 10x1 pipette |

#### **SLIM SG Pipettes**

|               |   |
|---------------|---|
| SG-SLIMX30-10 | ICSI SLIM Pipette with spike, ID 4-5µm, angle 30°, pack: 10x1 pipette |
| SG-SLIMX35-10 | ICSI SLIM Pipette with spike, ID 4-5µm, angle 35°, pack: 10x1 pipette |

#### **Holding SG Pipettes**

|            |  |
|------------|--|
| SG-H30a-10 | Holding Pipette, OD 90-110 µm, ID 15-20 µm, angle 30°, pack: 10x1 pipette  |
| SG-H35a-10 | Holding Pipette, OD 90-110 µm, ID 15-20 µm, angle 35°, pack: 10x1 pipette  |
| SG-H30b-10 | Holding Pipette, OD 110-130 µm, ID 20-23 µm, angle 30°, pack: 10x1 pipette |
| SG-H35b-10 | Holding Pipette, OD 110-130 µm, ID 20-23 µm, angle 35°, pack: 10x1 pipette |

#### **Biopsy SG Pipettes**

|            |   |
|------------|---|
| SG-B30B-10 | Blastomere Biopsy, Day 3, ID 35-42µm, angle 30°, pack: 10x1 pipette |
| SG-B35B-10 | Blastomere Biopsy, Day 3, ID 35-42µm, angle 35°, pack: 10x1 pipette |

|             |  |
|-------------|--|
| SG-B30Ta-10 | Blastocyst Biopsy, Day 5-6, Trophoectoderm, ID 25-30µm, angle 30°, pack: 10x1 pipette            |
| SG-B35Ta-10 | Blastocyst Biopsy, Day 5-6, Trophoectoderm, ID 25-30µm, angle 35°, pack: 10x1 pipette            |
| SG-B30Tb-10 | Blastocyst Biopsy, Day 5-6, Trophoectoderm, Universal, ID 30-35µm, angle 30°, pack: 10x1 pipette |
| SG-B35Tb-10 | Blastocyst Biopsy, Day 5-6, Trophoectoderm, Universal, ID 30-35µm, angle 35°, pack: 10x1 pipette |

|            |  |
|------------|--|
| SG-B30P-10 | Polar Body Biopsy, ID 10-15µm, angle 30°, pack: 10x1 pipette |
| SG-B35P-10 | Polar Body Biopsy, ID 10-15µm, angle 35°, pack: 10x1 pipette |

#### **Piezo Drill ICSI Pipettes**

|             |  |
|-------------|--|
| SG-XPi25-10 | Piezo Drill ICSI Pipette, ID 6 µm, angle 25° |
|-------------|--|