

SYNGA Microtools for IVF

CE certified, Gamma sterilized, MEA tested, nontoxic, safe for living cells.

Revision: 5-2026

I. Micromanipulation Pipettes

For use with a micromanipulator in IVF. Sterile, MEA, non-toxic, safe for living cells. Individually sealed pipettes, 10 pcs in a box.

ICSI Injection SG Pipettes

SG-X30S-10	ICSI Pipette with spike, OD 5-7µm, ID 4-5 µm, angle 30°
SG-X35S-10	ICSI Pipette with spike, OD 5-7µm, ID 4-5 µm, angle 35°
SG-X30WS-10	ICSI Pipette without spike, OD 5-7µm, ID 4-5 µm, angle 30°
SG-X35WS-10	ICSI Pipette without spike, OD 5-7µm, ID 4-5 µm, angle 35°

TESE SG Pipettes

SG-X30T-10	ICSI Pipette TESE, no spike, OD 7-9µm, ID 5-7µm, angle 30°
SG-X35T-10	ICSI Pipette TESE, no spike, OD 7-9µm, ID 5-7µm, angle 35°

SLIM SG Pipettes

SG-SLIMX30-10	ICSI SLIM Pipette with spike, OD 5-7µm, ID 4-5µm, angle 30°
SG-SLIMX35-10	ICSI SLIM Pipette with spike, OD 5-7µm, ID 4-5µm, angle 35°

Holding SG Pipettes

SG-H30a-10	Holding Pipette, OD 90-110 µm, ID 15-20 µm, angle 30°
SG-H35a-10	Holding Pipette, OD 90-110 µm, ID 15-20 µm, angle 35°
SG-H30b-10	Holding Pipette, OD 110-130 µm, ID 20-23 µm, angle 30°
SG-H35b-10	Holding Pipette, OD 110-130 µm, ID 20-23 µm, angle 35°

Biopsy SG Pipettes

SG-B30B-10	Blastomere Biopsy Pipette, Day 3, OD 40-47µm, ID 35-42µm, angle 30°
SG-B35B-10	Blastomere Biopsy Pipette, Day 3, OD 40-47µm, ID 35-42µm, angle 35°

SG-B30Ta-10	Blastocyst Biopsy Pipette, Day 5-6, Trophoctoderm, OD 27-33µm, ID 25-30µm, angle 30°
SG-B35Ta-10	Blastocyst Biopsy Pipette, Day 5-6, Trophoctoderm, OD 27-33µm, ID 25-30µm, angle 35°
SG-B30Tb-10	Blastocyst Biopsy Pipette, Day 5-6, Trophoctoderm, Universal, OD 33-38µm, ID 30-35µm, angle 30°
SG-B35Tb-10	Blastocyst Biopsy Pipette, Day 5-6, Trophoctoderm, Universal, OD 33-38µm, ID 30-35µm, angle 35°

SG-B30P-10	Polar Body Biopsy Pipette, OD 13-18µm, ID 10-15µm, angle 30°
SG-B35P-10	Polar Body Biopsy Pipette, OD 13-18µm, ID 10-15µm, angle 35°

Piezo Drill ICSI Pipettes

SG-X4Pi25-10	Piezo Drill ICSI ID 4	Piezo Drill ICSI Pipette, ID 4 µm, angle 25°
SG-X45Pi25-10	Piezo Drill ICSI ID 4,5	Piezo Drill ICSI Pipette, ID 4,5 µm, angle 25°
SG-X5Pi25-10	Piezo Drill ICSI ID 5	Piezo Drill ICSI Pipette, ID 5 µm, angle 25°
SG-X6Pi25-10	Piezo Drill ICSI ID 6	Piezo Drill ICSI Pipette, ID 6 µm, angle 25°
SG-X4Pi35-10	Piezo Drill ICSI ID 4	Piezo Drill ICSI Pipette, ID 4 µm, angle 35°
SG-X45Pi35-10	Piezo Drill ICSI ID 4,5	Piezo Drill ICSI Pipette, ID 4,5 µm, angle 35°
SG-X5Pi35-10	Piezo Drill ICSI ID 5	Piezo Drill ICSI Pipette, ID 5 µm, angle 35°
SG-X6Pi35-10	Piezo Drill ICSI ID 6	Piezo Drill ICSI Pipette, ID 6 µm, angle 35°

Assisted Hatching Needles

SG-AH30-10	Assisted Hatching Needle, angle 30°. Zona (ZP) opening, blastocoel collapsing (before vitrification)
SG-AH35-10	Assisted Hatching Needle, angle 35°. Zona (ZP) opening, blastocoel collapsing (before vitrification)